xiAPI.NET

xiAPI.NET

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Offline Processing

Documentation

xiAPI.NET stands for XIMEA Application Programming Interface for Dot Net - **Microsoft Visual C#**

It is a new common interface for all XIMEA cameras and represents a simplified version of generic **xiAPI**.

Architecture

API a is the software interface between the camera, operating system driver and application.

xiAPI.NET run-time is based on xiCam class implemented in *xiAPI.NET.dll* calling the *xiapi32.dll* on 32bit Windows operating systems and *xiAPI.NETX64.dll* calling the *xiapi64.dll* on 64bit Windows operating systems.

Installation

To add xiAPI.NET support:

- install **XIMEA Windows Software Package**
- When choosing components make sure: API and Examples are checked.
- Use one of the available versions of the .NET SDK for your application, supported versions are:
 - Microsoft .NET Framework 4.7.2
 - .NET Core 2.2
 - .NET Core 3.1
- We recommend to use Microsoft Visual Studio 2013 or a later version for application development.

Sample

.NET Sample code - See this section with various .NET examples.

After installation of the XIMEA Windows Software Package - the sample project is installed in the directory XIMEA/Examples/xiAPI.NET.

Here follows a sample code that captures ten images from the camera and saves them as .bmp files.

Exposure time is set 2 ms and gain is set to 5 dB.

```
// Sample for XIMEA API.NET

using System;
using System.Collections.Generic;
using System.Text;
using System.Collections;
using System.Drawing;
using xiApi.NET;
namespace xiApi.NET_example
{
class Program
{
static void Main(string[] args)
{
xiCam myCam = new xiCam();

try
{
// Initialize first camera
myCam.OpenDevice(0);

// Set device exposure to 2 milliseconds
int exposure_us = 2000;
myCam.SetParam(PRM.EXPOSURE, exposure_us);

// Set device gain to 5 decibels
float gain_db = 5;
myCam.SetParam(PRM.GAIN, gain_db);

// Set image output format to monochrome 8 bit
myCam.SetParam(PRM.IMAGE_DATA_FORMAT, IMG_FORMAT.MONO8);

//Start acquisition
myCam.StartAcquisition();
}
```

```

// Capture images
Bitmap myImage;
int timeout = 1000;
for (int i = 0; i < 10; i++)
{
    myCam.GetImage(out myImage, timeout );
    string fName = string.Format("image{0}.bmp", i);
    myImage.Save(fName);
}
// Stop acquisition
myCam.StopAcquisition();
}

catch (System.ApplicationException appExc)
{
    // Show handled error
    Console.WriteLine(appExc.Message);
    System.Console.ReadLine();
    myCam.CloseDevice();
}

finally
{
    myCam.CloseDevice();
}
} // end of Program
} // end of namespace xiApi.NET_example

```

Interface

The core of xiAPI.NET are the following functions, which allow the control of most of the camera functionality.

<code>// open interface</code>
<code>xiCam.OpenDevice(int DevId);</code>
<code>// set parameter</code>
<code>myCam.SetParam(PRM paramName, int paramVal);</code>
<code>myCam.SetParam(PRM paramName, float paramVal);</code>
<code>// get parameter</code>
<code>myCam.GetParam(PRM paramName, out int paramVal);</code>
<code>myCam.GetParam(PRM paramName, out float paramVal);</code>
<code>myCam.GetParam(PRM paramName, out string paramVal);</code>
<code>// get next image from buffer</code>
<code>myCam.GetImage(out myImage);</code>
<code>// close interface</code>
<code>xiCam.CloseDevice();</code>

Offline Processing

xiAPI.NET allows to process of already captured and stored images using Offline Processing.

Please read more at [xiAPI.NET Offline Processing](#).

Documentation

[xiAPI.NET Manual](#)

[Accessing xiAPI.NET pixel data in .NET Bitmap](#)