XIMEA CamTool

XIMEA CamTool

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Plugin development

CamTool description

- * The **XIMEA CamTool** is a cross-platform viewer application that is capable of setting basic camera parameters, controlling certain features, displaying the video stream and saving snapshots and sequences.
- * It is free to download and use. You do not need any activation.
- * This viewer is supported by all XIMEA camera families.
- * It runs on Windows, Linux, macOS operating systems offering a substantial imaging toolset that can be further extended with custom modules using a plugin infrastructure. **CamTool uses Qt under LGPL** for the UI and xiAPI for the camera control.
- * The CamTool's camera settings menu resembles the parameter set of the **xiAPI**.
- * It is part of **XIMEA Windows Software Package**.

- * Further, the program can be extended with custom functionality using its plugin infrastructure.
- * The default modules are already included and additional ones can be added using the plugin infrastructure.
- * CamTool supports the Processing chain feature (top right in the picture) where the user can see the pipeline of processes from the camera to view, and split it into multiple views.

Functions

The tool offers multiple functions like:

- * to see live image from multiple XIMEA cameras connected
- * control the camera parameters
- * store of camera image and video ([How to record video/image sequence inCamTool](#))
- * analyze the image properties
- * histogram and line profile
- * image averaging, image flip/mirror
- * software trigger timer, save/load camera and program settings
- * **LUT (Look up table)**
- * **Lua scripting**
- * **CUDA GPU Processing capabilities**

The modular plugin structure of the utility allows to add more and more features continuously.

CamTool supports the Processing chain feature (top right in the picture) where the user can see the pipeline of processes from the camera to view, and split it into multiple views.

CamTool description

The CamTool is a cross-platform viewer application that is capable of setting basic camera parameters, controlling certain features, displaying the video stream and saving snapshots and sequences.

Further, the program can be extended with custom functionality using its plugin infrastructure.

The default modules are already included and additional ones can be added using the plugin infrastructure.

CamTool is based on **Qt** for the UI and **xiAPI** for the camera control.


The CamTool's camera settings menu resembles the parameter set of the xiAPI.

Usage

- * Start CamTool with the executable can be found in C:\XIMEA\xiCamTool\.
- * Upon starting, CamTool will enumerate all connected cameras. If a camera is connected or disconnected during operation, a pop-up will appear which shows the name of the respective camera.
- * All available cameras are queued in the top bar menu, showing their names, acquisition controls and further information.
- * By default, only basic camera parameters are displayed in the camera settings menu. Further controls can be activated by switching the view mode to "Expert" or "Guru": Main Menu-> Edit->Options-> Visibility level
- * Each camera's image view can be moved and scaled using the left mouse button or scroll wheel.
- * **xiB** and **xiT** cameras have additional lens control function for the EF-mount: aperture and focus. The aperture can be set directly by selecting the respective F-value using the aperture slider. The focus requires movement value for the stepper motor to be defined first. Positive numbers will direct the movement to infinity. Negative numbers will direct the movement to macro. Then the specified move can be applied with the "Focus Move" button.
- * CamTool saves each camera's parameters and most user interface settings when the program is closed. The saved values are restored on the next start.

Adjusting the default settings

- * There are two ways how to reset cameras to the default setting in CamTool:
 - go to *Edit > Options* and in the dialog click the checkbox next to *Start with default camera parameter values next time* and click OK. After this quit CamTool and start it again.
 - or click on the small arrow next to the camera info icon and select *Restart Camera with Default Settings*


- * By default Auto-Exposure/Auto-Gain (AEAG) is enabled. Disable this feature using the “AEAG” checkbox to fine-tune the exposure and gain settings

Bandwidth Limitation

- * By default the camera uses an automatic bandwidth calculation which measures the available interface bandwidth when the device is started. After the measurement, a bandwidth limit is set using a 10-20% (this depends on the camera model) safety margin.
- * You can control the bandwidth limitation manually by changing the “Limit Bandwidth” value in the performance section. In order to evaluate the system limits, please monitor the ratio between triggered, transported and processed frames using the Tools->Camera Framerate view. If the processed frame rate is lower than the transported or triggered frame rate, it is required to reduce the bandwidth limit value.
- * The bandwidth limitation is especially helpful when setting up multi-camera systems to distribute the available bandwidth among all devices.


Multiple cameras

- * CamTool allows to connect and stream from multiple cameras and offers a few features for simplification of their control:
 - When more cameras are connected, labels with their names are displayed in the top panel.
When a label of a specific camera is clicked, its setting panel is raised to the top.
 - Cameras can be sorted alphabetically according to their names.
To activate it, check *Sort cameras alphabetically* in *Edit > Options > Cameras* :

Recording

- * CamTool allows to store the sequence of images in various formats. There are available multiple formats:
 - xiSeq, which is XIMEA proprietary format
 - **H.264**, AVI, which are standard formats
- * We test the playback of recorded data using the CamTool and standard formats also with **VLC**. Other decoders might work, but XIMEA does not guarantee full functionality in all camera modes.
- * Store speed is determined by the availability of CPU resources and the speed of the target drive. If you want to create a higher frame rate record, you can use the LoopRecord function that stores the image data first into the computer RAM. Later the result can be stored to other media and/or compressed by CPU.
See more details in the article: [Recording of Videos and Image Sequences in CamTool](#)

Recording of Videos and Image Sequences in CamTool

The CamTool enables to capture and store videos or image sequences in three Recording modes:

Record video to file - for recording directly to the drive


Record loop - for recording to memory or to drive in defined cycle length

Time-lapse - for recording image sequences with low to extremely low frame rates directly to the drive

All recording modes use the same file formats for storing videos or image sequences to the drive.

Supported recording file formats

CamTool supports saving to some of the most commonly used video formats and also in a sequence of individual frames.

Video formats that are supported with respective settings:

AVI file (*.avi)

Codec - OpenH264 or Uncompressed video can be selected

Quality - from 1 to 10, where 10 is the best image quality

H.264 file (*.h264)

Quality - from 1 to 10, where 10 is the best image quality

Note: We test the playback of recorded data using the CamTool and also with VLC. Other decoders might also work, but the full functionality in all camera modes is not guaranteed.

Image sequences that are supported with respective settings:

Numbered images sequence (*.xiseq) - image sequence is saved in separate folder beside *.xiseq file, which is a reference to this sequence. Image sequence can be stored in following file formats: *.bmp, *.jpg, *.png, *.gif, *.fits, *.tif, *.dng


TIFF multipage (*.tif)

Compression method - LZW, packbits or Uncompressed
enable/disable storing 16-bit data
enable/disable scale the color bit depth to a full 16-bit range

Record video to file

This mode allows to save video directly to the selected drive.

Shortcut for recording video to file is **CTRL + R**.

Pros

unlimited length of the record

Cons

recording frame rate is limited by hardware (writing speed of used drive) and by software (compression of stored data)


Record Loop

This recording mode allows to store captured images to memory or to drives. It allows to record from several cameras simultaneously (since V4.17.49).

When loop record comes to an end of the loop size, the recording process continues to store the next image by replacing the oldest recorded image.

This ensures that the last N images (frames) will be preserved after recording is stopped.


Shortcut for recording loop is **CTRL + L**.

It opens a dialog window:

Recording can be enabled for one or more cameras (checked in 1st column). Settings (storage, acq. buffer size...) are applied to the highlighted camera. More than one camera can be highlighted by holding the Ctrl button and clicking on unselected. The settings are common for all selected. Storage is the location of unprocessed transfer data. The recorded sequence can be saved afterward in the desired format. Only local drives are supported as storage.

When an error or skipped frame appears, the value is highlighted by red in the table.

Loop recording to memory

Pros

- highest possible frame rate
- unlimited duration of the recording

Cons

- limited number of recorded frames by the size of free memory

Loop recording to drive

Pros

- much longer recordings against loop recording to memory
- unlimited duration of the recording

Cons


- frame rate limited by writing speed of the drive

Size of the loop recording

The size of the loop recording determines how many recent frames will be stored. It can be set in a number of **frames**.

The maximum Loop size depends on image size and available memory or drive's space.

[Acq. buffer size defines the size of the circular buffer which contains image data from the sensor.](#)

Recording


Recording can be started with **START** button.

Recording can be stopped with **STOP** button. After stopping the recording process, all recorded frames are ready to be played in CamTool, or to be saved.

Saving of video or image sequence


For saving image sequence or video, press **Save as** button under the main menu of CamTool.

Shortcut for "Saving as" is **CTRL + S**.

Record Time-lapse

This recording mode allows to specify the interval between two consecutive frames which are to be stored to drive. Shortcut for recording time-lapse is **CTRL + P**.

Before Time-lapse is started, the number of images to be acquired has to be set. After acquisition reaches this image count, the acquisition is stopped.

Timing of Time-lapse

Two timing modes are available:

Free run - frame rate is calculated with actual settings of the camera

Interval - frame rate is forced - the interval between two consecutive images is set.

This interval has a span from milliseconds to hours

Output from Time-lapse

Two options are available:

Movie - image sequence can be stored in any supported video format

Custom sequence - image sequence can be stored in any supported image format

All stored image sequences or video files can have appended time stamp with format YYYYMMDD_HHMMSS_SSS. It is possible to append Coordinated Universal Time (UTC) instead of local time.

Performance

Tips for the best performance

For achieving the best performance during recording meet the following recommendations.

General tips

use xiCOP to diagnose potential performance issues

System Bandwidth

For the fastest performance of recording videos, the used system should have the highest possible bandwidth.

It can be influenced by:


type of used cable

slot, where the cable is connected (USB cameras)

number of connected cameras to one USB controller (USB cameras)

Measuring loop record writing speed

Measuring loop record writing speed is started with **Test write speed** button.

After pressing **STOP** button, the measured value is displayed in table.

Record Enable	Camera	Write speed [MB/s]	Transferred images
<input checked="" type="checkbox"/>	MC124CG-SY-UB (CUCAU1718009)	9284	0

Type	Path
RAM	

These measured results show approximate average writing speed to the selected storage.

The actual writing speed might differ.

Limitations

Recording speed is limited by settings of the camera:

- exposure time
- image format
- image height (ROI)
- bandwidth limit
- device output bit depth

File format AVI is limited to size 511GB.

Image Modes

- * RAW8: 1 byte (8 bit) per pixel transport and buffer format
RAW16: 10 or 12 bit transport in packed or unpacked format, 2 bytes (16 bit) buffer format with padding
 - Unpacked: every pixel is transported as 2 bytes, which results in half the transfer speed, compared to RAW8
 - Packed: every pixel is combined with the following pixel to the next bigger number of bytes, e.g. 24x 12bit in 192bit. If the packed mode is used, during the call of `xiGetImage`, the image data is unpacked to 2 bytes per pixel. Since this processing step is executed on the CPU and requires certain resources, depending on the image size, the overall acquisition frame rate will be lower than the measured transport speed.
- * RGB24: RAW8 image data is transported and then converted into RGB (demosaicing) during the call of `xiGetImage`. The result is 3 channels R, G; B with 8 bits per channel.
- * RGB32: Similar to RGB24 but with an additional 4th channel, representing the alpha value. This mode results in larger images, but usually faster processing.
- * Mono8/Mono16: On mono cameras equivalent to RAW8/RAW16. On color cameras, the result is a single channel grey image, after the demosaicing (interpolation) has been applied.
- * TRANSPORT: This mode provides access to the transported image data, without further processing, such as unpacking or demosaicing. One possible application is direct storage of the packaged 12bit data for later 'offline' processing.

Annotations and Scale bar

CamTool enables adding simple text annotations to the image. It also supports a scale bar. In order to work with annotations, you need to show the tool using menu Tools - Annotations. You might find more details [here](#).


Temperature plugin

While working with cameras, developers and customers find it useful to have a graph representing known camera temperatures.

Here are the functions and features of the XIMEA plugin:

- * The plugin reads temperature values from all temperature sensors located in various parts of the camera and displays them in the graph window.
- * The graph shows the change in temperature during the time since the "Start measure" button was pushed. The "Stop measure" button ends the measurement.
- * The window also informs whether the cooling is active.
- * The sampling period of the temperature is adjustable. For example: every 5 seconds.
- * The plugin allows to export the formatted data to the clipboard.
- * The temperature is recorded simultaneously from all connected cameras, but only one dataset is presented in the graph - it comes from the camera selected in the Combo box at the top.

- * Time can be reset by the Reset All button. This also clears all the data previously gathered so please be careful.
- * Autoscale - Automatic horizontal and vertical scale. Possible is zoom and drag in different directions.

Additional Controls:

The temperature control can be switched on or off in the general control section on the right.

The target temperature value can be also set in the general control section.

Note: Plugin included in XIMEA Software package since V4.19.24

Formats supported

Various types of files that can be selected as output during Save as:

- BMP**
- JPEG**
- PNG**
- GIF**
- DNG**
- TIFF**

Download, Install and Run

- * Windows: XIMEA CamTool is part of [XIMEA Software Package](#). [Download](#)
 - Run by clicking on XIMEA CamTool icon  created on Desktop
- * macOS - XIMEA CamTool is part of [XIMEA macOS Software Package](#). [Download](#)
 - Run by clicking on XIMEA CamTool icon  created in Applications
- * Linux - XIMEA CamTool is part of [XIMEA Linux Software Package](#)
 - Run by starting shell command `/opt/XIMEA/bin/xiCamTool`

Minimal System Requirements

Windows

Microsoft Windows 10,11 (64-bit)

CPU: Intel Core i5

Available RAM: 2GB for view operations

Available RAM: 3GB for recording with compression

Available Hard drive space: 400MB

macOS


macOS version 11.3 (Big Sur)

CPU: Intel Core i5 or Apple M1

Available RAM: 3GB

Available Hard drive space: 300MB

The following error refers to minimal system requirements:

Plugin development

It is possible to implement a custom plugin for the XIMEA CamTool. More information can be found [here](#).

Feedback can be sent to info@ximea.com